

Basic of Proximate Analysis and New Soxtec 8000 and Hydrocap 8000

Presented by Yaovalak Suksong
Sithiporn Associates Co.,Ltd.



Basic of Proximate Analysis and New Soxtec 8000 and Hydrocap 8000

- ▶ **Part 1 : Basic of Proximate Analysis**
- ▶ Part 2 : New Soxtec 8000 and Hydrocap 8000

▶ 2

SAC 2013

Purpose of Proximate Analysis

Estimation and determination of how much of the major food components, which are Moisture, CHO, Lipids, Proteins, Ash, Crude Fibre, exist in a given food. The proximate analyses therefore are:

1. Moisture
2. Crude Fat
3. Crude Protein
4. CHO and Crude Fibre

Total carbohydrate = 100-[moisture + crude fat + crude protein + ash].

▶ 3

SAC 2013

- ✱ **Sample Preparation**
- ✱ **Crude Protein – Kjeldahl Method**
- ✱ **Fibre**
 - ✱ **Crude Fibre**
 - ✱ **Detergent Fibre**
 - ✱ **Dietary Fibre**
- ✱ **Fat**
 - ✱ **Crude Fat**
 - ✱ **Total Fat**

▶ 4

SAC 2013

Sample Preparation



Feature	Cemotec™	Cyclotec™	Knifetec™	Homogenizers
Sample Type	Dry samples prior to moisture analysis, up to 15% Moisture and 10% Fat, such as: Grains, Cereals Seeds Beans Dry granular foods Pelleted feed Fertilizer	Dry samples prior to wet chemistry or IR analysis, up to 15% Moisture and 10% Fat, such as: Grains, Cereals Seeds Petfood pellets Hay, straw (dried and cut to 2-3cm) Silage (dried and cut to 2-3 cm) Pelleted feed Leaves Lime, Coal, Tobacco	High-moisture, high-fat and fibrous samples, such as: Grains, Cereals Beans Oilseeds Nuts Peas - wet and dry Maize (Corn) - wet and dry Pelleted feed and Pet food (3-4 mm) Pellets up to 6 mm - with lid for pellets Meat products Vegetables and Fruit	High-moisture, high-fat and fibrous samples, such as: Meat and meat products Whole prepared meals Forage, hay, straw and silage (2097 recommended) Fish* and fish products Vegetables and fruit Chemical and pharmaceutical formulations <i>*May require descaling and or skinning</i>
Sample Size	Up to 14 mm Ø	Up to 10 mm, large inlet up to 40 mm	Maximum 100 ml (50 - 150 g)	(2094): 0,1 - 1,5 kg (2097): 0,1 - 2.5 kg
Grinding Principle	Two discs, one rotating one stationary	Impeller, abrasive ring, and screen	Rotor blade	Various rotor blade

▶ 5

**Crude Protein
Kjeldahl**



▶ 6

SAC 2013

Kjeldahl Method



✚ **Johan Kjeldahl** → Chemist in Copenhagen, Denmark

✚ **Kjeldahl Procedure**


- ✚ **Digestion** => the reaction between organic compounds and sulfuric acid produced **ammonium salts**.
- ✚ **Distillation** => Ammonium salts reacted with strong alkali. the **ammonia** produced in this step was distilled and dissolved in a standardized solution of hydrochloric acid or sulfuric acid.
- ✚ **Titration** => the solution was back titrated with Sodium hydroxide to indirectly measure nitrogen.

✚ During the 1880s, *Kjeldahl used potassium sulfate to raise the boiling point of the acid and mercury as a catalyst to speed the decomposition.* For the back titration process of the released ammonia, he used boric acid buffer solution.

▶ 7 http://en.wikipedia.org/wiki/Johan_Kjeldahl

SAC 2013

ISO 5983-1 Animal Feeding stuffs-Determination of nitrogen content and calculation of crude protein content Part 1 :




- **Wight sample** → 0.5-2g
- **Digestion** →
 - Add 15 g $K_2SO_4 + 0.9-1.2g CuSO_4$
 - Add 25 ml H_2SO_4 for 1g sample and 6-12 ml for each additional gram of sample
 - Digest 2h
- **Distillation** →
 - Add water 250-350 ml
 - 100 ml 33% NaOH
 - Receiver solution 25 ml
 - ✚ 0.05M, 0.125M H_2SO_4
 - ✚ 4% H_3BO_3 solution
 - Steam distill until ≥ 150 ml of distillate
- **Titration** →
 - ✚ 0.05M, 0.125M H_2SO_4 → 0.1M, 0.25M NaOH
 - ✚ 4% H_3BO_3 solution → 0.05M, 0.125 H_2SO_4
- **Calculation** →
 - Nitrogen Content
 - % Protein

▶ 8

SAC 2013

The Kjeldahl Procedure



▶ 9

SAC 2013

Sample preparation of Nitrogen analysis

- **Solid Sample**

- The particle size should be equal to, or less than, 1 mm

- **Semi Solid Samples**

- Mortar

- **Liquid Samples**

- The sample contain precipitates filtration or sedimentation and the filtrate or supernatant liquid taken for analysis
- Shaking or stirring

- **Sample Weight**

- Analytical balance accurate to 0.1 mg
- The actual weight of sample required is dependent on Nitrogen content
 - ▶ Homogenous samples (excluding water) 0.1 - 1.0 g
 - ▶ Non homogenous samples 1.0 - 3.0 g or more
 - ▶ Water samples (dependant on N content) 1.0 - 100 ml




▶ 10

SAC 2013

The Digestion Step




Breakage of nitrogen bonds and reduction to Ammonium ions (NH_4^+)



▶ 11

SAC 2013

The Digestion System




▶ 12

SAC 2013





The Distillation Step

Separation of ammonia (NH₃)



Digestion mixture
plus
Water 80ml (6:1)
40% NaOH 50-60ml (4:1)



▶ 13

SAC 2013

FOSS-Kjeltec™ 8000 series



Kjeltec 8420



Kjeltec 8460

▶ 14

SAC 2013

Kjeltec Safety – SAFE Patented


Delay

- ▶ Dilution
- ▶ Alkali addition
- ▶ Delay 2 - 20 seconds to allow for acid/alkali reaction
- ▶ Distillation



SAFE (Steam Addition for Equilibration)

- ▶ Dilution
- ▶ Distillation starts - mixing the acid and water to ensure dilution takes place
- ▶ Gradual alkali addition during distillation slowly raises the pH minimizing reaction



New Global standard
EN ISO 5983-2
2005

▶ 15

SAC 2013

Titration



- Titration with standard acid determines the amount of ammonia and therefore nitrogen in the sample.
- Protein can be calculated from a known nitrogen content.



Initial color



Plus NH₃



Intermediate



Final end point

Receiver Solution
Boric acid + indicators
(Methyl Red &
Bromocresol Green)

Titration with titrant acid
(HCl or H₂SO₄)

▶ 16

SAC 2013

Titration & Calculation



Automatic Titrator



Digital Burette



Standard Glassware Burette

$$\% \text{ N} = \frac{(T - B) \times N \times 14.007}{\text{wt of sample (mg)}} \times 100$$

$$\% \text{ Protein} = \% \text{ N} \times \text{factor}$$

T = Volume of titrant used for Sample (ml)

B = Volume of titrant used for Blank (ml)

N = Normality of titrant (to 4 decimal places)

14.007 = Molecular weight of Nitrogen

▶ 17

SAC 2013

Quality Control – Recovery



Digestion

- Standard substance (Tryptophan, Acetanilide)
- Certified Reference sample
- Internal Reference sample

Distillation

- (NH₄)₂SO₄ (purity ≥ 99.5%)
- (NH₄)₂Fe(SO₄)₂·6H₂O

▶ 18

SAC 2013



☹ **Low/high results**

☹ **Poor Repeatability**

▶ 19

SAC 2013

Trouble Shooting

Sample Preparation

- ☑ Balance
- ☑ Sample size
- ☑ Homogeneous

Digestion

Incomplete digestion

- Volume of acid (H_2SO_4), salt (K_2SO_4), catalyst (Cu, Hg, Se, Ti, Cu/Ti)
- Temperature of block
- Digestion Time
- Exhaust head, aspiration rate
- Vacuum level

▶ 20

SAC 2013

Trouble Shooting

Digestion

▶ Foaming Sample

- Reduce sample size
- Boiling Rods
- Anti-foam agents → n-Octanol, H_2O_2
- Ramp the temperature



▶ Salt cake formation

- Add a small amount (20-30 ml) of water to the sample when cool enough to handle
- Reheat the mixture in a block to “melt” the cake
- Distill the sample with SAFE.....

▶ 21

SAC 2013

Trouble Shooting

Distillation and Titration

- ▶ Alkali volume
- ▶ Distillation recovery
 - Dry ammonium sulfate before use → 103-105 °C
 - Purity ≥ 99.5% Titrant concentration
- ▶ Burette, air bubbles
- ▶ Adsorption of nitrogen in distilling unit (carry over effect), dirty splash head
- ▶ An uncontrolled blank contribution, nitrogen in the chemicals/reagents used

- ▶ **Acid/Base Mixture Problem**
 - Distill the sample with SAfE

▶ 22

SAC 2013

Application



DIGESTION SYSTEMS

- ▶ Kjeldahl digestion
- ▶ Chemical Oxygen Demand (C.O.D.)
- ▶ Trace and Heavy metal preparation
- ▶ Sample preparation for CFA, FIA etc.



DISTILLATION SYSTEMS

- ▶ Kjeldahl Nitrogen/Protein
- ▶ Ammonia
- ▶ Nitrate
- ▶ Nitrite
- ▶ Sulphite
- ▶ Phenols
- ▶ Cyanides
- ▶ Volatile Acidity
- ▶ etc


▶

SAC 2013

Conclusion

Kjeldahl Method

- ❖ Digestion
- ❖ Distillation
- ❖ Titration



Nitrogen /Protein



Trouble shooting

▶ 24

SAC 2013 July 18, 2013

Fibre



▶ 25

SAC 2013

Synopsis



- What is fibre?
- Fibre analysis
 - ⊗ Crude fibre
 - ⊗ Detergent fibre
 - ⊗ Dietary fibre
- Trouble shooting



▶ 26

SAC 2013

What is fibre?



Fibre contains cellulose, hemicellulose and lignin + other substances that are difficult or impossible to digest.



▶ 27

SAC 2013

Fibre Analysis



- ▶ **Crude Fibre - Feed trading**
- ▶ **Neutral Detergent Fibre - Optimisation of Energy intake**

➡ Fibertec 2010/M6



- ▶ **Dietary Fibre – Food labelling**

➡ Fibertec E



▶ 28


SAC 2013

Sample Preparation – Crude Fibre

- ▶ A particle size < 1 mm.
- ▶ Homogeneous
- ▶ **Fat content**
 - Fat content > 10% → defatted prior to analysis.
 - Fat content is between 1% and 10%
 - ❖ defatting is recommended (ISO5498)
 - ❖ should be performed (AOAC, AACC, AOCS)
 - Fat content is less than 1% defatting is not necessary.
- ▶ Defat with **Petroleum Ether**
- ▶ **High Carbonate**
 - Carbonate (CO_3^{2-}) content is above 1 % carbonate
 - 0.5 M HCl
- ▶ **Sample Weight**
 - An analytical balance accurate to 0.1 mg

▶ 29


SAC 2013



▶ 30

SAC 2013

Crude & Detergent fibres



▶ 31

SAC 2013

Steps in Crude Fibre analysis



AOAC 4.6.01 - Filter Method

- ▶ Weigh into beaker
- ▶ Defatting if > 1 % fat
- ▶ **TRANSFER** to beaker
- ▶ Boil with 1.25 % acid
- ▶ Filtration and wash
- ▶ **TRANSFER** to beaker
- ▶ Boil with 1.25 % alkali
- ▶ Filtration and wash
- ▶ Wash with alcohol
- ▶ **TRANSFER** to ashing disc
- ▶ Ashing

Fibertec procedure

- ▶ Weigh into crucible
- ▶ Defatting - Cold Extraction Unit
- ▶ Crucible to Fibertec
- ▶ Boil with 1.25 % acid
- ▶ Filtration and wash
- ▶ Boil with 1.25 % alkali
- ▶ Filtration and wash
- ▶ Wash with acetone
- ▶ Ashing

▶ 32

SAC 2013

Transfer of Samples Potential Problems



- ⊗ Manual work → Takes time
- ⊗ Loss of sample → Low result
- ⊗ Exposure to reagents → Hazard




▶ 33

SAC 2013

Crude Fibre – Weende Method



DEFINITION
The residue of plant cells after extraction by Acid and Alkaline Hydrolysis




▶ 34

SAC 2013

Detergent Fibre – Van Soest Method



DEFINITION
The residue of plant cells after fractionation using detergent solutions.




NDS = Neutral Detergent Solution
ADS = Acid Detergent Solution

▶ 35

SAC 2013

Fibre Calculation



● = CH
▲ = Fibre
■ = Protein

$$\% \text{ fibre} = \frac{W2 - W3}{W1} \times 100$$

▶ 36

SAC 2013

Residue Content



<p>❏ Crude fibre</p>	<p>Cellulose 50-80% Hemicellulose ~20% Lignin 10-50%</p>
<p>❏ Detergent fibre</p>	
<p>Neutral Detergent Fibre (NDF)</p>	<p>Cellulose 100% Hemicellulose 100% Lignin 100%</p>
<p>Acid Detergent Fibre (ADF)</p>	<p>Cellulose 100% Lignin 100%</p>
<p>Acid Detergent Lignin (ADL)</p>	<p>Lignin 100%</p>

▶ 37

SAC 2013

Conclusion – Crude and Detergent Fibre



Method	Principle	Residue content
Crude Fibre (Weende method)	Crude fibre is defined to be the residue after sequential treatment with hot H ₂ SO ₄ (conc. 1,25 %) and hot NaOH (conc. 1,25 %)	Cellulose, 50 - 80% Hemicellulose, approx 20% Lignin, 10 - 50%
NDF (acc. to van Soest)	Neutral Detergent Fibre is defined to be the residue after treatment with a neutral detergent solution (Sodium lauryl sulphate and EDTA)	Cellulose, 100% Hemicellulose, 100% Lignin, 100%
ADF (acc. to van Soest)	Acid Detergent Fibre is defined to be the residue after treatment with an acid detergent solution (Cetyl trimethylammonium bromide in Sulphuric acid solution)	Cellulose, 100% Lignin, 100%
ADL (acc. to van Soest)	Acid Detergent Lignin is defined to be the residue after initial treatment by the ADF method followed by removal of the cellulose fraction through extraction using 72% H ₂ SO ₄	Lignin, 100%

▶ 38


SAC 2013

Dietary Fibre



DEFINITION

Dietary fibre is the remnants of plant components that are resistant to hydrolysis by the digestive enzymes of man



▶ 39

SAC 2013

Fibertec E



Fibertec E



▶ 40

SAC 2013

Sample preparation – Dietary Fibre



■ Solid Sample

- particle size 0.3- 0.5 mm



■ Semi solid

- ▶ Homogenize and dry samples prior to grinding.
- ▶ If sample cannot be heated, freeze-dry before milling. The sample

■ High fat → fat content > 10 % (recommended (ISO 5498))

- Defat with 25 ml **Petroleum Ether** for three times


■ High Carbohydrate → Sugar > 50%

- Extract three times each with 10 ml of **Ethanol**

▶ 41

SAC 2013

Determination of TDF, IDF & SDF



▶ 42

SAC 2013

Trouble Shooting



📌 Repeatability

📌 Low/high results

▶ 43

SAC 2013

Trouble Shooting



📌 Sample Preparation

- ☺ Homogeneous
- ☺ Crude Fiber
 - ☺ Size ≤ 1 mm.
 - ☺ Fat $> 1\%$, Carbonate $> 1\%$
- ☺ Dietary Fiber
 - ☺ Size 0.3-0.5 mm
 - ☺ Fat $> 10\%$, Sugar $> 50\%$

📌 Crucible

- ☺ Clean → Chromic acid, Hydrochloric acid solution (15-20%)

▶ 44

SAC 2013

Trouble Shooting



📌 Celite

- ☺ accuracy ± 0.1 mg

📌 Foaming

- ☺ Antifoam → n-octanol
- ☺ Turn down heater
- ☺ Check cooling water

📌 Wash sample

- ☺ Hot water

▶ 45

SAC 2013

Conclusion



- **Crude Fiber**
 - Fibertec M6, Fibertec 2010
- **Detergent Fiber**
 - Fibertec M6, Fibertec 2010
- **Dietary Fiber**
 - Fibertec E
- **Quality Control**
 - Trouble Shooting

▶ 46

SAC 2013

Solvent Extraction Fat



▶ 47

SAC 2013

Topic

- ▶ Methods for fat determination:
 - ▶ Direct Extraction → Crude Fat
 - ▶ Total Fat
- ▶ Quality Control
 - ▶ Trouble Shooting
 - ▶ Suggestion



▶ 48


SAC 2013

Principle – Fat Analysis




Crude Fat

- ▶ Crude fat is determined by a solvent extraction.
- ▶ non-polar organic solvents such as hexanes, petroleum ether.



Total Fat


- ▶ Total fat determination includes a preparatory acid hydrolysis step and a solvent extraction
- ▶ The fat that is bound to other non-solvent solubles as e.g. proteins are separated in hydrolysis step
- ▶ Hydrolysis makes chemically or mechanically bound fats accessible to solvent extraction



▶ 49

SAC 2013

Solvent Extraction: The Soxtec idea - to improve the Soxhlet idea



- | | |
|----------------|--|
| Faster | -An automatic four step extraction
-Loading/unloading of six samples simultaneously |
| Cheaper | -Less solvent use
-Solvent recovery
-Batch operation |
| Safer | -Explosion proof design
-Less exposure to solvent vapours |

▶ 50

SAC 2013

Sample Preparation



- **Solid Sample**
 - The particle size should be equal to, or less than, 1 mm
- **Semi Solid Samples**
 - Homogenizing == > Mortar
- **Liquid Samples**
 - Samples containing particles, like wastewater, are filtered to collect the parts that will be extracted
 - In samples where it is a suspension, handling including adsorption on inert material such as Celite 566 could be the



Fat Content	Recommended Sample Weight
0-10%	1.5 – 2g ± 0.1 mg
10 – 20%	1 – 1.5g ± 0.1 mg
>20%	1g ± 0.1 mg

▶ 51

SAC 2013





Crude Fat and Total Fat Analysis



Crude Fat	Total Fat
<ul style="list-style-type: none"> ▶ Sample Preparation ▶ Solvent Extraction - Soxtec 	<ul style="list-style-type: none"> ▶ Sample Preparation ▶ Hydrolysis – 3-6N HCl
	 <ul style="list-style-type: none"> ▶ Solvent Extraction - Soxtec

Soxtec Principle



			
<p>Boiling Step</p> <ul style="list-style-type: none"> ▪ Sample is immersed in boiling solvent ▪ Provides rapid extraction of soluble materials ▪ For most application, the boiling step is 15 to 25 minutes 	<p>Rinsing Step</p> <ul style="list-style-type: none"> ▪ Sample is raised out of the boiling solvent ▪ Condensed solvent drips through sample and rinses out residuals ▪ Usually 30 to 40 minutes 	<p>Solvent Recovery Step</p> <ul style="list-style-type: none"> ▪ Condensed solvent is collected in the collection vessel ▪ This step concentrates the extracted material in the extraction cup and saves solvent for reuse ▪ Usually 10 minutes 	<p>Pre-drying Step</p> <ul style="list-style-type: none"> ▪ Shortened drying time through pre-drying in instrument

Soxtec System





NEW



Soxtec 8000

Hydrotec 8000



Soxtec 2055



Soxtec 2043



Soxcap 2047



Soxtec 2045



☹ **Low/high results**

☹ **Repeatability**

▶ 55

SAC 2013

Trouble Shooting



☺ **Check sample size**

➡ Size \leq 1 mm.

☺ **Homogeneous**

➡ Celite

☺ **Sample weight**

➡ Accuracy \pm 0.1 mg

➡ Negative weight

☺ **Cup temperature**

➡ Adjust level between cup and heater

☺ **Drop rate**

➡ 3-5 drops/s

▶ 56

SAC 2013

Recommendation



Before start.....

Check

- Instrument
 - Soxtec
 - Cooling
- Apparatus
 - Cup
 - Thimble

Sample Preparation

- After Hydrolysis \rightarrow Dry sample

Solvent

- AR Grade

▶ 57

SAC 2013

Application



Soxtec

Sample

- Soil, wastewater and sludge.
- Plastics, petrochemicals, paper, textiles and a very wide range of other industrial matrix

Parameters

- Extractable Matter – Material soluble in a given solvent or range of solvent

Soxtec + Hydrotec/Soxcap

Sample

- Raw Materials, Intermediates and Finished Products in Food, Animal Feed and Petfood

Parameters

- Total Fat (Free and Bound Fat)

▶ 58

SAC 2013



Part 2 : New Soxtec 8000 and Hydrocap 8000



SAC 2013
Sithiporn Scientific & Technology Conference

New SOXTEC 8000 & HYDROTEC 8000
Solvent Extraction Solution

Presented by Miss Patcharapom Paruang
Sithiporn associates co.,Ltd

FOSS

**Soxtec 2050,2055,2043,2045
&
Soxcap 2047**



**The new Soxtec 8000
&
Hydrotec 8000**

Dedicated Analytical Solutions

Product Intro
Reinventing Fat Analysis – 3 new products



Hydrotec™ 8000



Hydrocap




Soxtec™ 8000

Hydrocap	Patented unique Hydrocap single use filter system
Hydrotec 8000	New fully automated hydrolysis (12 pos.)
Soxtec 8000	New fully automated expandable extraction unit (6+6 pos.)

Dedicated Analytical Solutions

First fully automated Soxhlet system

- Patented Hydrocap system for integrated total fat analysis
- Solvent dial for enclosed solvent dosage
- Intelligent water cooling, pivoting hotplates



HydroCap contains sample all the way through



Dedicated Analytical Solutions


- Save operator time by avoiding filter transfer in total fat
- Flexible extraction according to demand (84 sample/day)
- Simplified operation reduces risk of errors and associated costs



- *Avoid sample handling*
- *Save water*
- *Save solvent*
- *Improve throughput*

Dedicated Analytical Solutions

- Automation reduces risk of contact with chemicals and solvents
- Technology simplifies operation and reduces the risk of accidents
- Safe system allows broad range of solvents



- Safe External control

- No acid handling
- Automated acid collection

- Speed dial for closed solvent dosage
- Effective cooling to avoid evaporation
- Solvent recovery with leakage sensors

Dedicated Analytical Solutions

- Easy to install (optional customer installation)
- Individual hotplates per cup/sample (run less than 6 sample&no cup adjustment after docking)
- No magnets or thimble adapters made of metal
- Control Unit manage one/two Extraction units (6 or 12 positions system)
- Improved capacity:
 - 7 batches allows for 42 analyses/day or 84/day
 - one batch overnight
- Improved user Interface with big display and graphics
- Solvent leakage sensor and solvent sensor in tank to prevent start of analysis if tank is not emptied
- Improved handling of solvent drainage (flask design and drainage system)
- Backward compatibility – Soxcap 2047 or Hydrolysis Unit 1047



Dedicated Analytical Solutions

Hydrotec 8000

- Automated hydrolysis
- Acid filling, boiling, cooling, rinsing, draining
- Closed hydrolysis system
- eliminates emission of acidic corrosive fumes
- Batch handling – with up to 12 samples at a time
- 9 user defined programs



Hydrocap

- Unique Hydrocap single use filter system
- Same sample holder all the way
- from weighing to extraction
- Hydrocap fits into any Soxhlet system



Dedicated Analytical Solutions



Dedicated Analytical Solutions